



Sil-Fixx Sterile Transparent Dressing



Conlett Sil-Fixx Sterile Transparent Dressing uses silicone adhesive technology.

Sil-Fixx Sterile is a Semi-Permeable dressing with a micro-porous transparent film. This provides a high moisture vapour transmission rate, protecting against fluid contaminants and bacteria.

Features of the Sil-Fixx Sterile Transparent Dressing:

- Transparent Semi-Permeable Dressing
- Micro-Porous transparent film with high moisture vapour transmission rate
- Protects against fluid contaminants and bacteria.
- Easy application
- Conforms to the skin easily, very little roll up incidence
- Ideal for, paper skin, neonatal, elderly, dialysis, chemotherapy, & radiography patients.
- Pain free and stress free removal

Indications:

Intended for use on IV, donor, incision, after removal of IV catheter, vaccine sites, minor burns, ulcers, skin tears, lacerations and minor wounds.

Contraindications:

It is not recommended that the products are to be applied over deep wound cavities, second or third degree burns or where the wound shows evidence of infection.

SIZES AVAILABLE

When ordering please quote the following

DESCRIPTION	DIMENSIONS	ORDER#	QTY/BOX
Sterile Transparent Dressing	4cm x 6cm	SFCM0406	50
Sterile Transparent Dressing	4cm x 10cm	SFCM0410	50
Sterile Transparent Dressing	6cm x 7cm	SFCM0607	50

(Minimum order quantity of 1 box)

Supramed

U25/17 Cairns St
Loganholme
QLD 4129
Australia

Tel: (+61) 7 3209 9957

Email: contact@supramed.com.au