



# Sil-Fixx Sterile Transparent Dressing



Conlett Sil-Fixx Sterile Transparent Dressing uses silicone adhesive technology.

Sil-Fixx Sterile is a Semi-Permeable dressing with a micro-porous transparent film. This provides a high moisture vapour transmission rate, protecting against fluid contaminants and bacteria.

## Features of the Sil-Fixx Sterile Transparent Dressing:

- Transparent Semi-Permeable Dressing
- Micro-Porous transparent film with high moisture vapour transmission rate
- Protects against fluid contaminants and bacteria.
- Easy application
- Conforms to the skin easily, very little roll up incidence
- Ideal for, paper skin, neonatal, elderly, dialysis, chemotherapy, & radiography patients.
- Pain free and stress free removal

## Indications:

Intended for use on IV, donor, incision, after removal of IV catheter, vaccine sites, minor burns, ulcers, skin tears, lacerations and minor wounds.

## Contraindications:

It is not recommended that the products are to be applied over deep wound cavities, second or third degree burns or where the wound shows evidence of infection.

## SIZES AVAILABLE

When ordering please quote the following

DESCRIPTION	DIMENSIONS	ORDER#	QTY/BOX
Sterile Transparent Dressing	4cm x 6cm	SFCM0406	50
Sterile Transparent Dressing	4cm x 10cm	SFCM0410	50
Sterile Transparent Dressing	6cm x 7cm	SFCM0607	50
Sterile Transparent Dressing	10cm x 12cm	SFCM1012	50

(Minimum order quantity of 1 box)

## Conlett

U25/17 Cairns St  
Loganholme  
QLD 4129  
Australia

Tel: (+61) 7 3209 9957  
Email: [contact@conlett.com](mailto:contact@conlett.com)